

SAHPRA/2024/DESTRUCTION OF MEDICINES, MEDICAL DEVICES AND IVDS/RFB009

Questions and Answers

Questions	Answers
<p>1. May you please assist with the following, what type of packaging are you referring to please specify. e.g 5-liter, 10 litre, or 20 litres? For collection</p>	<ul style="list-style-type: none"> • The medicines or products that will be collected are usually in bulk containers already from the companies. Some medicines or products will have been collected from operations with or without law enforcement authorities and these will be sealed in forensic bags (SAPS or SAHPRA ones). The quantity of the consignment is usually large (in Kilograms and can even reach 1000kg). The type of packaging material needed should be able to prevent spillage or breakages of the products from the storage side (Point of origin) to the destruction site. Huge bags may be used, or the cartons may be wrapped for easy carriage in the transportation vehicle. Big boxes (cartons) may be used also.